# POLYMAILLE® C

- FR** Prothèses vasculaires en polyester tricoté
- EN** Knitted polyester vascular prostheses
- DE** Gefäßprothesen aus gestricktem Polyester
- ES** Prótesis vasculares de poliéster tejido
- IT** Protesi vascolari in poliesteri intrecciato
- PT** Próteses vasculares em poliéster tricotado
- NL** Vaatprothesen uit gebreide polyester
- PL** Protezy naczyńiowe z poliestru trykotowego
- DA** Vaskulære proteser af strikket polyester
- SV** Av polyester stickade kärlproteser
- NO** Vaskulære proteser i strikket polyester
- FI** Polyesteritrikoiset verisuoniproteesit
- EL** Αγγειακά μοσχεύματα από πλεκτό πολυεστέρα
- TR** Örülü polyester vasküler protezler
- RO** Proteze vasculare din poliester tricostat
- HU** Kötött poliészter érprotézisek
- RU** Протезы сосудов из вязаного полиэстера
- CS** Cévní protézy z pleteného polyesteru
- AR** توصف الوصلات الوعائية من البوليستر المحبوك
- BG** Съдови протези от плетен полиестер

**STERILE R**

FR : stérilisé par irradiation  
EN : sterilized by irradiation  
DE : sterilisiert durch Bestrahlung  
ES : esterilizado por irradiación  
IT : sterilizzato con irradiazione  
PT : sterilise pela irradiação  
NL : gesteriliseerd per uitstraling  
PL : sterylizacji metodą napromieniowania  
DA : steriliseret ved bestråling  
SV : steriliserats med strålning  
NO : sterilisert ved bestråling  
FI : teriloitu säteilyty  
EL : ασποραιομένο από την ακτινοβολία  
TR : ışınlama ile steril  
RO : sterilizat prin iradiere  
HU : sterilizálni besugárzás  
RU : простерилизовано облучением  
CS : sterilizovaný zářením  
AR : تعقيمها من اشعاع  
BG : стерилизирани чрез облъчване

FR : Lire la notice  
EN : See instruction for use  
DE : Anleitung lesen  
ES : Consulte las instrucciones de uso  
IT : Leggere attentamente le istruzioni per l'uso  
PT : Consulte as instruções de utilização  
NL : Lees de handleiding  
PL : Zobacz instrukcję obsługi  
DA : Se brugsanvisning  
SV : Se bruksanvisningen  
NO : Les bruksanvisningen  
FI : Lue käyttöohje  
EL : Βλ. οδηγίες χρήσης  
TR : Kullanma kilavuzuna bakınız  
RU : См. инструкцию по использованию  
UK : Дивіться інструкцію з використання  
CS : Viz návod k použití  
RO : Consultați instrucțiunile de utilizare  
HU : Lásd a használati útmutatót  
AR : الإطلاع على نشرة التعليمات  
BG : Вижте указанията

FR : Usage unique,  
EN : Single use,  
DE : Zum einmaligen Gebrauch  
ES : De un solo uso  
IT : Monouso  
PT : Utilização única  
NL : Eenmalig gebruik  
PL : Do jednorazowego użytku  
DA : Engangsbrug  
SV : Engångsbruk  
NO : Engangsbruk  
FI : Kertakäyttö  
EL : Μιας χρήσης  
TR : Tek kullanımlık  
RU : Одноразового применения  
UK : Одноразового використання  
CS : Na jedno použití  
RO : De unică folosință  
HU : Egyszeri használatra  
AR : يُستخدم مرة واحدة  
BG : За еднократна употреба

FR : Conserver à l'abri de la lumière du soleil  
EN : Keep away from sunlight  
DE : Unter Lichtabschluss aufbewahren  
ES : Conservar en un lugar protegido de la luz solar  
IT : Conservare al riparo dalla luce del sole  
PT : Manter abrigado da luz solar  
NL : Beschermen tegen het zonlicht  
PL : Przechowywać z daleka od słońca  
DA : Beskyttes mod sollys  
SV : Skyddas från solljus  
NO : Lagres beskyttet fra sollys  
FI : Säilytettävä suojassa auringonvalolta  
EL : Να φυλάσσεται μακριά από το ηλιακό φως  
TR : Güneş ışığından uzak tutunuz  
RU : Хранить в месте, защищенном от солнечных лучей  
UK : Зберігати в місці, яке захищено від сонячних променів  
CS : Zabraňte přístupu slunečního světla  
RO : A se păstra la adăpost de lumina soarelui  
HU : Napfénytől távol tartandó  
AR : يُحفظ بعيداً عن أشعة الشمس  
BG : Да се пази от пряка слънчева светлина

FR : Ne pas utiliser si l'emballage est endommagé  
EN : Do not use if the packaging is damaged  
DE : Nicht verwenden, wenn die Verpackung beschädigt ist  
ES : No utilizar si el embalaje está deteriorado  
IT : Non utilizzare se la confezione è danneggiata  
PT : Não utilizar se a embalagem se encontrar danificada  
NL : Niet gebruiken als de verpakking beschadigd is  
PL : Nie używać jeśli opakowanie jest uszkodzone  
DA : Må ikke anvendes, hvis emballagen er beskadiget  
SV : Får inte användas om förpackningen är skadad  
NO : Skal ikke brukes hvis emballasjen er skadet  
FI : Ei saa käyttää, jos pakkaus on vioittunut  
EL : Να μην χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά  
TR : Ambalaj hasar götmüşse kullanılmayınız  
RU : Не использовать в случае повреждения упаковки  
UK : Не використовувати у разі пошкодження упаковки  
CS : Nepoužívejte, pokud je porušen obal  
RO : A nu se utiliza dacă ambalajul este deteriorat  
HU : Ne használja fel, ha a csomagolás sérült  
AR : لا يُستخدم في حالة وجود أي تلف بالعبوة  
BG : Да не се употребява, ако е нарушена целостта на опаковката.

FR : Craint l'humidité  
EN : Protect from humidity  
DE : Vor Feuchtigkeit schützen  
ES : Proteger de la humedad  
IT : Teme l'umidità  
PT : Manter a humidade  
NL : Kan niet tegen vocht  
PL : Trzymać z daleka od wilgoci  
DA : Beskyttes mod fugt  
SV : Skyddas från fukt  
NO : Beskyttes mot fuktighet  
FI : Suojattava kosteudelta  
EL : Να προστατεύεται από την υγρασία  
TR : Nemi koruyunuz  
RU : Поддерживайте алажность  
UK : Підтримуйте вологість  
CS : Udržujte vlhkost  
RO : A se proteja de umezeală  
HU : Szárraz helyen tartandó  
AR : يُحفظ بعيداً عن الرطوبة  
BG : Да се пази от влага

**LOT**

FR : Numéro de lot  
 EN : Batch code  
 DE : Chargennummer  
 ES : Número de lote  
 IT : Numero di lotto  
 PT : Número de lote  
 NL : Lotnummer  
 PL : Numer partii produktu  
 DA : Batch-kode  
 SV : Lotnummer  
 NO : Lotnummer  
 FI : Eränumero  
 EL : Κωδικός παρτίδας  
 TR : Parti kodu  
 RU : Номер партии  
 UK : Номер партії  
 CS : Císlo šarže  
 RO : Număr lot  
 HU : Adagszám  
 AR : رقم التسلسل  
 BG : Партида номер

**SN**

FR : Numéro de série  
 EN : Serial number  
 DE : Seriennummer  
 ES : Número de serie  
 IT : Numero di serie  
 PT : Número de série  
 NL : Serienummer  
 PL : Numer serii  
 DA : Serienummer  
 SV : Serienummer  
 NO : Serienummer  
 FI : Sarjanumero  
 EL : Αριθμός σειράς  
 TR : Seri numarası  
 RU : Номер серии  
 UK : Номер серії  
 CS : Sériové číslo  
 RO : Număr serie  
 HU : Sorozatszám  
 AR : رقم التسلسل  
 BG : Серийен номер

FR : À utiliser avant la date  
 EN : Use by  
 DE : Verfallsdatum  
 ES : Utilizar antes de  
 IT : Utilizzare entro  
 PT : Utilizar antes de  
 NL : Te gebruiken vóór  
 PL : Do użycia przed  
 DA : Anvendes inden  
 SV : Används före  
 NO : Brukes før  
 FI : Käytettävä ennen  
 EL : Χρήση έως  
 TR : Son kullanma tarihi  
 RU : Использовать до  
 UK : Використовувати до  
 CS : Použitelné do  
 RO : A se utiliza înainte de  
 HU : Minőségét megőrzi:  
 AR : تاريخ انتهاء الصلاحية  
 BG : Да се използва преди

**REF**

FR : Numéro de référence  
 EN : Reference number  
 DE : Bestellnummer  
 ES : Número de referencia  
 IT : Numero di catalogo  
 PT : Referência de catálogo  
 NL : Catalogusnummer  
 PL : Numer katalogowy  
 DA : Referencenummer  
 SV : Referensnummer  
 NO : Henvisningsnummer  
 FI : Viitenumero  
 EL : Αριθμός αναφοράς  
 TR : Referans numarası  
 RU : Номер по каталогу  
 UK : Номер за каталогом  
 CS : Referenční číslo  
 RO : Număr de referință  
 HU : Hivatkozási szám  
 AR : كود المنتج  
 BG : Каталоген номер

FR : Emballage recyclable  
 EN : Recyclable packaging  
 DE : Verpackung wiederverwertbar  
 ES : Embalaje reciclable  
 IT : Confezione riciclabile  
 PT : Embalagem reciclável  
 NL : Recycleerbare verpakking  
 PL : Opakowanie do użytku wtórnego  
 DA : Emballagen kan genanvendes  
 SV : Förpackningen kan återvinnas  
 NO : Resirkulerbar emballasje  
 FI : Kierrätettävä pakkaus  
 EL : Ανακυκλώσιμη συσκευασία  
 TR : Geri dönüştürümlü ambalaj  
 RU : Упаковка подлежит вторичной переработке  
 UK : Упаковка підлягає вторинному переробленню  
 CS : Recyklovatelný obal  
 RO : Ambalaj reciclabil  
 HU : Újrahasznosítható csomagolás  
 AR : عبوة قابلة لإعادة التدوير  
 BG : Опаковка за рециклиране

FR : Ne pas restériliser  
 EN : Do not sterilise  
 DE : Nicht resterilisieren  
 ES : No volver a esterilizar  
 IT : Non risterilizzare  
 PT : Não voltar a esterilizar  
 NL : Niet hersterilisieren  
 DA : Ma ikke resteriliseres  
 PL : Nie poddawać ponownie sterylizacji  
 SV : Skall inte återanvändas  
 NO : Skal ikke brukes på nytt  
 FI : Ei saa uudestaansteriloida  
 EL : Να μην επαναποστεριώνεται  
 TR : Yeniden sterilize etmeyin  
 RO : Nu reesterilizați  
 HU : Ne sterilizálja újra  
 RU : Не стерилизовать  
 CS : Znovu nesterilizujte  
 AR : لا تقم بإعادة التعقيم  
 BG : Не извършвайте повторна стерилизация

FR : Limite de température  
EN : Temperature limit  
DE : Temperaturbegrenzung  
ES : Limite de temperatura  
IT : Limiti di temperatura  
PT : Limite de temperatura  
NL : Temperatuurbepanking  
DA : Temperaturgrænse  
PL : Ograniczenie temperatury  
SV : Temperaturgräns  
NO : Temperaturgrense  
FI : Lämpötilaraja  
EL : Όριο θερμοκρασίας  
TR : Sıcaklık limiti  
RO : Limită de temperatură  
HU : Hőmérséklet határérték  
RU : Предельная температура  
CS : Teplotní limit  
AR : حدود الحرارة  
BG : Температурен лимит

FR : Fabricant  
EN : Manufacturer  
DE : Hersteller  
ES : Fabricante  
IT : Fabricante  
PT : Fabricante  
NL : Fabrikant  
DA : Producent  
PL : Producent  
SV : Tillverkare  
NO : Produsent  
FI : Valmistaja  
EL : Κατασκευαστής  
TR : İmalatçı  
RO : Producător  
HU : Gyártó  
RU : Производитель  
CS : Výrobce  
AR : المصانع  
BG : Производител

FR : Consulter les précautions d'emploi  
EN : Consult instructions for use  
DE : Bitte die Vorsichtsmaßnahmen beachten  
ES : Véanse precauciones de uso  
IT : Consultare le istruzioni per l'uso  
PT : Consultar as instruções de utilização  
NL : Raadpleeg de voorzorgen bij gebruik  
DA : Se brugsanvisningen  
PL : Zapoznać się ze środkami ostrożności  
SV : Se försiktighetsåtgärder  
NO : Se forholdsregler ved bruk  
FI : Tutustu varotoimenpiteisiin  
EL : Συμβουλευτείτε τις οδηγίες χρήσης  
TR : Kullanım talimatlarına bkz.  
RO : Consultați instrucțiunile de utilizare  
HU : Olvassa el a használati utasításokat  
RU : Ознакомьтесь с мерами предосторожности  
CS : Přečtěte si pokyny pro používání  
AR : ضرورة الرجوع لاحتياطات الاستعمال  
BG : Вижте указанията за употреба

## INSTRUCTIONS ET INDICATIONS

### I/ DESCRIPTION

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## *INSTRUCTIONS AND INDICATIONS*

### **I/ DESCRIPTION**

POLYMAILLE®C vascular prostheses are constructed from knitted polyester. The models available are impregnated with collagen of bovine origin.

Collagen impregnation of the POLYMAILLE®C vascular prostheses ensures water permeability of less than 5 ml/cm<sup>2</sup>/min at a pressure of 120 mm of mercury (Hg) and allows direct implantation of the prosthesis without the need for pre-clotting, thereby reducing operating time and blood loss.

### **II/ INDICATIONS**

POLYMAILLE®C vascular prostheses are indicated for replacement or bypass of arteries presenting aneurysm or obliterative arterial disease. Their indication is restricted to abdominal and peripheral surgery not crossing the knee flexion crease.

### **III/ CONTRAINDICATIONS**

POLYMAILLE®C vascular prostheses are contraindicated for surgery of the thoracic aorta, coronary bypass, arteriovenous access, extra-corporeal circulation and crossing the knee flexion crease.

POLYMAILLE®C vascular prostheses are contraindicated in patients with an immune reaction or hypersensitivity, known or suspected, to collagen of bovine origin.

### **IV/ POSSIBLE COMPLICATIONS**

The most commonly encountered complications following implantation of polyester vascular prostheses include, but are not restricted to: infections, periprosthetic collections, thromboses, pseudoaneurysms, periprosthetic haematomas, blood loss, revascularisation oedema.

A thrombectomy on a prosthesis can lead to risks of damaging the collagen coating.

In the event of direct puncture of the prosthesis during endovascular procedures, there is a risk of damaging the prosthesis.

### **V/ WARNINGS AND PRECAUTIONS**

The warnings listed below must be strictly observed. Failure to observe them may result in significant blood loss, a risk of loss of a limb, or death of the patient.

☛ **POLYMAILLE®C vascular prostheses are for single use only. DO NOT REUSE.** Any part of the prosthesis remaining after the procedure must be destroyed. POLYMAILLE®C prostheses should never be saved for reimplantation in another patient or resteralised by any method whatsoever.

After use or removal, the prostheses must be considered biologically hazardous and non-functional and must be handled and discarded in accordance with medical practices and the local, national and federal laws in force.

## INSTRUCTIONS AND INDICATIONS

- **POLYMAILLE®C vascular prostheses are sterilised by irradiation. DO NOT RESTERILIZE.**
- **DO NOT USE IF THE PACKAGE IS DAMAGED.** 
- **DO NOT USE POLYMAILLE®C VASCULAR PROSTHESES AFTER THEIR EXPIRY DATE.**
- Store protected from light, cold, heat ( $10^{\circ}\text{C} < T < 35^{\circ}\text{C}$ ) and humidity ( $\text{RH} < 80\%$ ).
- **DO NOT PRE-CLOT THE POLYMAILLE®C VASCULAR PROSTHESES** to avoid any risk of infection. They are impregnated with collagen.
- Prostheses contaminated with blood and/or other body fluids prior to surgery must not be used.
- Atraumatic clamps should be used or the jaws of the clamps should be enclosed in soft protective covers. Tighten as little as possible. There is a risk that impaired collagen coating may cause a leak.
- **AVOID EXCESSIVE MANIPULATION OF THE PROSTHESIS IN THE EVENT OF PRECLOTTING AND DURING ITS IMPLANTATION.** There is a risk of damaging the collagen coating which may cause a leak.
- **AVOID EXCESSIVE STRETCHING AND AXIAL TORSIONS (twisting)**, which may result in kinking.

- Suturing should be performed using tapered needles only; sharp-edged needles may damage the prosthesis.

## VI/ INSTRUCTIONS FOR USE

- Vascular prostheses should only be inserted in the operating theatre by surgeons trained in vascular surgery techniques.
- On removal from their cardboard packaging, the POLYMAILLE®C prostheses appear as follows: the peel-open aluminium foil pouch contains a double pack (an outer blister within which a sterile second peel-open blister contains the prosthesis). Carefully open the first (external) blister and entrust the internal sterile blister to persons in sterile dress. This should be opened by the surgical team. Use new sterile gloves to handle the prosthesis.

## *INSTRUCTIONS AND INDICATIONS*

### **ATTENTION: DO NOT OPEN THE ALUMINIUM POUCH ONTO THE STERILE DRAPE.**

This pouch is simply a protective package and its contents are not sterile.

- The diameter of the prosthesis must correspond to the diameter of the host vessel. It is important to make allowance for the patient's weight, the operating position and the degree of movement likely to occur at the implantation site when determining the correct length of the prosthesis.

### **VII/ TRACEABILITY**

Self-adhesive labels and a patient card are provided in each package to allow the required follow-up of implanted prostheses. One label must be pasted onto the patient's medical record and another onto the patient card.

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