1997

 **PEROUSE**
MEDICAL

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VIII/ INFORMATION PATIENT

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I/ DESCRIPTION

POLYTHESE® vascular prostheses are woven polyester vascular grafts. Three models are available:

- POLYTHESE® IC and ICT, impregnated with collagen of bovine origin;
- POLYTHESE® IC-3GL, impregnated with collagen of bovine origin;

Collagen impregnation of POLYTHESE® IC, IC-3GL and ICT vascular prostheses ensures a water permeability of less than 10 ml/cm²/min at 120 mmHg and allows direct implantation of the prosthesis without the need for pre-clotting, thereby reducing operating time and blood loss.

POLYTHESE® IC-3GL prostheses have 3 equidistant longitudinal guide lines.
FORMALDEHYDE/GLUTARALDEHYDE-FREE

II/ INDICATIONS

POLYTHESE® vascular prostheses are indicated for replacement or bypass of arteries depicting aneurysm or arterial disease.

Their indication is restricted to thoracic (POLYTHESE® ICT), abdominal and peripheral surgery not involving the crossing of the knee joint.

The POLYTHESE® IC-3GL prostheses are especially indicated in aortic root and ascending thoracic aorta surgery.

III/ CONTRAINDICATIONS

POLYTHESE® vascular prostheses are not indicated for coronary bypass or arteriovenous access.

POLYTHESE® IC, IC-3GL and ICT prostheses are contraindicated for patients with known hypersensitivity to collagen.

IV/ POSSIBLE COMPLICATIONS

The most commonly encountered complications following implantation of polyester vascular prostheses include, but are not restricted to: infection, periprosthetic fluid collection, thrombosis, pseudoaneurysms, periprosthetic haematoma, leaking, rupture, lymphatic complications, venous thrombosis, pulmonary embolism.

V/ WARNINGS AND PRECAUTIONS

The afore-mentioned warning messages are to be strictly respected; for the patient, their non-compliance can result in significant blood loss, disability of a limb or even death.

- ☛ **POLYTHESE® vascular prostheses are for single use only. DO NOT REUSE. POLYTHESE® prostheses must never be reimplanted. As they are contaminated with blood, a used vascular prosthesis constitutes a biological risk (infectious risk). In addition, PEROUSE MEDICAL cannot guarantee proper functioning of the device after explantation (water-tightness of the prosthesis; integrity of the collagen coating; weakening of the textile structure, particularly at the areas of anastomosis; level of burst strength; level of suture retention strength, etc.)**
- ☛ **POLYTHESE® vascular prostheses are sterilised by irradiation. DO NOT RESTERILIZE.**
- ☛ **DO NOT USE THE IMPLANT IF THE PACKAGE IS DAMAGED.**
- ☛ **DO NOT USE POLYTHESE® VASCULAR PROSTHESES AFTER THEIR EXPIRY DATE.**

- ☛ Store away from light, cold, heat (10°C < T < 35°C) and humidity (RH < 80%).
- ☛ **DO NOT PRE-CLOT POLYTHESE® IC, IC-3GL AND ICT PROSTHESES: they are impregnated with collagen.**
- ☛ Prostheses contaminated with blood and/or other body fluids prior to surgery must not be used.
- ☛ Do not use scissors to cut the prosthesis to the desired length.
- ☛ The use of a thermocutter can prevent fraying of the extremities of the prosthesis.
- ☛ Atraumatic clamps should be used or the jaws of the clamps should be enclosed in soft protective covers. Tighten as little as possible.
- ☛ Avoid excessive tension on the prosthesis during implantation.
- ☛ Avoid excessive length and axial torsions (twisting), which may result in kinking.
- ☛ Suturing should be performed using only tapered needles; sharp-edged needles may damage the prosthesis.

VI/ INSTRUCTIONS FOR USE

- ☛ Vascular prostheses should only be implanted by surgeons trained in vascular surgery techniques.
- ☛ On removal from their outer packaging, the prostheses are as follows: peel-open aluminium pouch containing the double packaging (external blister pack containing a second sterile blister pack holding the prosthesis).
WARNING: DO NOT OPEN THE ALUMINIUM POUCH ON THE STERILE DRAPE. It is simply a protective package and its contents are not sterile.
Peel open the aluminium pouch then carefully open the first blister back and hand the contents only to people wearing sterile clothing. The second blister pack is sterile and should be opened by the surgical team. New sterile gloves should be used to handle the prosthesis.
- ☛ The diameter of the prosthesis must correspond to the diameter of the patient's vascular lumen. The patient's weight must be considered the operating position and the degree of movement likely to occur at the implantation site when determining the correct length of the prosthesis to be implanted.

VII/ TRACEABILITY

Self-adhesive labels and a patient card are provided in each packaging to allow the required follow-up of implanted prostheses. One label must be stuck onto the patient's medical record and another onto the patient card.

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