


# Surpass™

SUPER-STIFF GUIDE WIRE



CROSSING BY BEST SUPPORT

AN ACCESSORY FOR PERIPHERAL INTERVENTION



Lifetech  
SCIENTIFIC

CE 0344



# Surpass™

## SUPER-STIFF GUIDE WIRE

The Surpass™ super-stiff guide wire is designed to provide maximum support for navigation of tortuous anatomy during peripheral intervention procedures.

### Easy to use


- A radiopaque gold inner coil allows precise tip visibility and location during wire-guide manipulation.

### Low friction surface

- Surpass™ providing maximum support and low friction surface with PTFE coating technology during complex diagnostic and interventional procedures.

### Soft and shapable tip

- Super long tip design stiff shaft and smooth transition at the junction to reduce transfer resistance and vascular damage.



Specification	Guidewire diameter (inch)	Guidewire length L1 (mm)	Guidewire inner spring length L2 (mm)
SSGW-35-260	0.035	2600	40

### Lifetech Scientific (Shenzhen) co.,Ltd.

Add: Cybio Electronic Building, Langshan 2nd Street, North Area of High-tech Park, Nanshan District, Shenzhen 518057, P. R. China

Tel: +86-755-86026250 Fax: +86-755-86026251

Email: marketing@lifetechmed.com

Web: www.lifetechmed.com

EU authorized representative:  
 LifeTech Scientific (Europe) Coöperatief U.A.  
 Oliemolenstraat 60, 6416CB Heerlen, the Netherlands