

## MOZEC™ NC TECHNICAL SPECIFICATIONS

<b>Catheter System</b>	: Rapid Exchange (Rx)
<b>Balloon Material</b>	: Nylon - Non-Compliant
<b>Available Balloon Lengths (mm)</b>	: 8, 10, 13, 15, 18, 23, 28, 30, 35, 38
<b>Available Balloon Diameters (mm)</b>	: 2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.00, 4.50
<b>Usable Catheter Length</b>	: 142 cm
<b>Max Guidewire Compatibility</b>	: 0.014" (0.36 mm)
<b>Guide Catheter Compatibility</b>	: 5F (min I.D. 0.056" / 1.42 mm)
<b>Distal Shaft Coating</b>	: Biocompatible, Hydrophilic & Lubricious coating from distal balloon base upto rapid exchange port
<b>Balloon RO Markers</b>	: 2 - Platinum / Iridium
<b>Tip Length</b>	: 3.5±0.5 mm
<b>Proximal Shaft Diameter</b>	: 1.98 F
<b>Distal Shaft Diameter</b>	: 2.7 F
<b>Markers on Proximal Shaft</b>	: Brachial and Femoral markers 90 cm and 100 cm from distal tip
<b>Nominal Pressure (NP)</b>	: 12 atm for all Diameters
<b>Rated Burst Pressure (RBP)</b>	: 20 atm for all Diameters

## MOZEC™ NC ORDERING INFORMATION

Dia / Length	8mm	10mm	13mm	15mm	18mm
2.00mm	MNC20008	MNC20010	MNC20013	MNC20015	MNC20018
2.25mm	MNC22508	MNC22510	MNC22513	MNC22515	MNC22518
2.50mm	MNC25008	MNC25010	MNC25013	MNC25015	MNC25018
2.75mm	MNC27508	MNC27510	MNC27513	MNC27515	MNC27518
3.00mm	MNC30008	MNC30010	MNC30013	MNC30015	MNC30018
3.50mm	MNC35008	MNC35010	MNC35013	MNC35015	MNC35018
4.00mm	MNC40008	MNC40010	MNC40013	MNC40015	MNC40018
4.50mm	MNC45008	MNC45010	MNC45013	MNC45015	MNC45018

Dia / Length	23mm	28mm	30mm	35mm	38mm
2.00mm	MNC20023	MNC20028	MNC20030	MNC20035	MNC20038
2.25mm	MNC22523	MNC22528	MNC22530	MNC22535	MNC22538
2.50mm	MNC25023	MNC25028	MNC25030	MNC25035	MNC25038
2.75mm	MNC27523	MNC27528	MNC27530	MNC27535	MNC27538
3.00mm	MNC30023	MNC30028	MNC30030	MNC30035	MNC30038
3.50mm	MNC35023	MNC35028	MNC35030	MNC35035	MNC35038
4.00mm	MNC40023	MNC40028	MNC40030	MNC40035	MNC40038
4.50mm	MNC45023	MNC45028	MNC45030	MNC45035	MNC45038

**CAUTION:**  
This product is intended for use by or under the direction of a Physician. Prior to use, it is important to read the package insert supplied with this device, thoroughly for Instructions for Use, Warnings and Potential Complications associated with use of this device. Tests performed by and data on file at Meril Life Sciences Pvt. Ltd. All drawings are artists representations only and not to be considered as an engineering drawing or photographs. Information for the use only in countries with applicable health authority product registrations.



**Meril Life Sciences Pvt. Ltd.**  
Survey No. 135/139, Bilakhia House,  
Muktanand Marg, Chala,  
Vapi 396 191.  
Gujarat, India.  
T +91 260 305 2100  
F +91 260 305 2125

**Meril, Inc.**  
116, Village Boulevard, Suite 200,  
Princeton, NJ 08540, USA.  
T +1 609 951 2287  
F +1 609 520 1702

**Meril South America**  
Doc Med LTDA  
Al. dos Tupiniquins,  
1079 - Cep: 04077-003 - Moema,  
Sao Paulo, Brazil.  
T +55 11 3624 5935  
F +55 11 3624 5936

**Meril GmbH.**  
Bornheimer Strasse 135-137,  
D-53119 Bonn,  
Germany.  
T +49 228 7100 4000  
F +49 228 7100 4001

**Meril Tibbi Cihazlar**  
İmalat Ve Ticaret A. S.  
6. Mimar sinan Mah.,  
Çavusbasi Cad. Özde Sok.  
Aydın Eksi İş Merkezi Kat:1  
Cekmekoy/Istanbul, Turkey.  
T +90 53 2272 5172

E askinfo@merillife.com  
W www.merillife.com

Mozec™nc is registered trademark  
of Meril Life Sciences Pvt. Ltd.

MNC/BROCHURE/003/MLS/20160902/IND

**MOZEC™ NC**  
Rx PTCA Balloon Dilatation Catheter

Meril

**MOZEC™ NC**  
Rx PTCA Balloon Dilatation Catheter

# MOZEC™ NC

Rx PTCA Balloon Dilatation Catheter

Handicraft for *Angiocraft*.

Experience controlled and confident expansion with MOZEC™ NC

Soft tip made up of medical grade PEBAX to facilitate advancement of the catheter to and through the stenosis

Low tip entry profile for easy cross and re-cross

Short abrupt balloon shoulders  
2 folds for Ø 2.00 mm  
3 folds for Ø 2.25 to 4.50 mm

Hump free fusion of distal tip inner tubing and balloon

Biocompatible, hydrophilic, lubricious coating from distal balloon base upto rapid exchange port

Nylon balloon material provides non-compliance at high pressure

Low crossing profile and highly flexible MOZEC™ NC is designed to acquire maximum pressure and meet clinical needs

Controlled balloon growth for flat compliance at high pressure

